

## Venturi valves and Eco mask kits

58% reduction in environmental impact<sup>1</sup>



## Venturi valves and kits now available in the Eco range

Introducing a comprehensive range of venturi valves available individually or in kits. These products are designed to meet the requirements of BS EN 13544-3: 2001 and colour coded for easy identification.

### Comfortable for the patient and the environment

Utilisation of the latest manufacturing technology has enabled us to combine two non-PVC materials in the same mask.

The material forming the body of the mask is clear and rigid enough to maintain the mask's shape. A second, softer material is utilised in the manufacture of the seal, which is in contact with the patient's face.

The use of these materials has resulted in a mask kit with an environmental impact score of 7.6 milli ecopoints when compared to the equivalent PVC mask kit, which has an environmental impact score of 18.1 milli ecopoints which represents a 58% reduction<sup>1,2</sup>.

### More products, more choice

Wherever you see the distinctive green logo you can be assured that we are working hard to lower the environmental impact of these products.



1. Adult Eco Mask with Venturi Valve environmental evaluation, IQR159 - 15630 N.Goodman, August 2008.  
2. 2006 SimaPro Version 6, Pre Consultants bv, Plotterweg 12, 3821 BB Amersfoort, The Netherlands

## Flow rates

Oxygen flow meter - litres per minute (L/min)

	2	3	4	5	6	7	8	9	10	11	12	15
24%	52	78	104	130	156	181	207	233	259	285	311	389
28%	22	34	45	56	67	78	90	101	112	123	135	168
31%	16	24	31	39	47	55	63	71	79	87	94	118
35%	11	17	23	28	34	39	45	51	56	62	68	84
40%	8	12	17	21	25	29	33	37	41	46	50	62
60%	4	6	8	10	12	14	16	18	20	22	24	30

**NOTE:** Figures in white indicate total gas flows which may be less than the patient's peak inspiratory flow.

Total flow rate to patient - litres per minute (L/min)

## Optimum flow rates

The flow printed on the valves ensures that the patient receives a total flow of approximately 40 L/min<sup>3</sup>. This is designed to meet the patients 'normal' peak flow requirements.


If the patient is hyperventilating then their peak flow requirements may exceed 40L/min<sup>3</sup>. The oxygen flow can be increased leading to a higher total flow (see table). The oxygen concentration will remain unchanged.

3. Clinical Practice in Respiratory Care. James B.Fink. P277


Try the product  
visit [www.intersurgical.com/ecorange](http://www.intersurgical.com/ecorange)

## Features and benefits

Using new non-PVC materials meant a new mask design was required, and this provided an ideal opportunity to make improvements to the traditional aerosol mask design.




## Ordering information



Code: 1107080

Box qty. 25


Adult Eco venturi kit



Code: 1040013

Box qty. 40

40% oxygen recovery kit



Code: 1024080

Box qty. 50

Adult Eco mask with 24% venturi valve



Code: 1028080

Box qty. 50

Adult Eco mask with 28% venturi valve



Code: 1031080

Box qty. 50

Adult Eco mask with 31% venturi valve



Code: 1035080

Box qty. 50

Adult Eco mask with 35% venturi valve



Code: 1040080

Box qty. 50


Adult Eco mask with 40% venturi valve



Code: 1060080

Box qty. 50

Adult Eco mask with 60% venturi valve



Code: 0024

Box qty. 50

24% oxygen venturi valve

Code: 0028

Box qty. 50

28% oxygen venturi valve

Code: 0031

Box qty. 50

31% oxygen venturi valve

Code: 0035

Box qty. 50

35% oxygen venturi valve

Code: 0040

Box qty. 50

40% oxygen venturi valve

Code: 0060

Box qty. 50

60% oxygen venturi valve