


Respi-Check™

The immediate and simple visual indicator of respiratory rate



The immediate and simple visual indicator of respiratory rate

Why measure respiratory rate?

Respiratory rate should be recorded in anyone requiring oxygen therapy, whether due to illness, trauma or recovery from sedation. Yet it is often inaccurately estimated and poorly recorded.

How do I measure Respiratory Rate?

The Respi-Check™ provides accurate real-time information about a patient's work of breathing. It comprises a visible ball contained within a clear plastic tube which has been integrated into a high concentration oxygen mask.


The flow of oxygen into the mask pushes the ball to the top of the tube. The slightest negative pressure generated by the patient's inspiratory effort pulls the ball to the bottom of the tube providing an instant visual notification of each breath taken.



Ordering information

The complete Respi-Check product range

Product	Code	Box Qty.
Adult high concentration mask Oxygen mask and tube	1202	20
Adult medium concentration mask Oxygen mask	1216	20
Paediatric high concentration Oxygen mask and tube	1292	20



IS10.7 • Issue 6 10.11