


## The route to humidification



A range of products to help you along the way



# The route to humidification - your choice

At Intersurgical we understand that each patient and clinical situation is different. Our aim is to provide a solution for all of your humidification requirements, active or passive, for both ventilated and spontaneously breathing patients.

## Ventilated

### Passive humidification

Passive Humidification requires a heat and moisture exchanger to be positioned at the patient connection of a basic two limb breathing system. This is designed to replicate the functions of the upper airway conserving the patient's own expired heat and moisture and returning these to the patient during inspiration.




### Active humidification

Active humidification requires a water bath humidifier, humidification chamber and either a heated wire breathing system or water traps within the system. This technique is also designed to replicate the functions of the upper airway by the addition of heat and moisture from the humidifier. This provides a higher level of humidity than an HME and should be selected depending on the clinical requirements of the patient.



## Spontaneous



## Why is Humidification needed

In normal respiration the upper airway helps to warm and humidify inspired air, and to retain the warmth and moisture contained in expired air. During inspiration even cold or dry air is typically heated to 37 °C and fully saturated, containing 44mg H<sub>2</sub>O per litre.

In mechanical ventilation the patient's upper airway may be bypassed by the introduction of a tracheal tube. As a result the patient's lungs may be confronted with dry inspired gas.


The drying and cooling effect is exacerbated by the presence of the tracheal tube, the normal process of re-absorption of heat and moisture by the upper airway during expiration is lost.

Prolonged exposure to dry ventilatory gases can lead to a number of problems as highlighted below.

### Prolonged exposure to dry ventilatory gases can lead to

- Localised inflammation of the trachea.
- A reduction in ciliary function
- Retention and thickening of secretions
- Lowering of patient temperature
- Reduction in Cardiopulmonary function
- Increased risk of tracheostomy tube occlusion

### - respiratory Epithelium adversely affected by heat moisture loss



## High routes are available

There are two options for patient humidification, **passive** or **active**.

**Passive humidification** conserves the patient's own heat and moisture whilst **active humidification** adds additional heat and moisture via a humidifier.

We have a full range of products for both options to suit your patients requirements.





# Passive Humidification

If your ventilated patient requires passive humidification then we have a wide range of basic breathing systems and heat and moisture exchangers.

In this brochure are a selection of options however, for the full range please refer to our website [www.intersurgical.com](http://www.intersurgical.com) or our product catalogue.

## Heat and Moisture Exchangers

Heat and moisture exchangers are designed to be used at the patient connection of a breathing system to prevent heat and moisture loss when the upper airway is by passed.

We can provide both HME only options and a range combined with filters, HMEFS. These provide the performance of a dedicated HME with the filtration efficiency of a breathing filter.

## Filta-Therm Plus and Clear-Therm HMEFs

Filtration efficiency	Moisture return at: VT500ml	Resistance to flow		Compressible volume	Weight	Minimum tidal volume
		30L/min	60L/min			
99.999	31.5mg H <sub>2</sub> O/L	1.3cm H <sub>2</sub> O	3.0cm H <sub>2</sub> O	66ml	44g	200ml
Filtration efficiency	Moisture return at: VT500ml	Resistance to flow		Compressible volume	Weight	Minimum tidal volume
		30L/min	60L/min			
99.99	32mg H <sub>2</sub> O/L	1.0cm H <sub>2</sub> O	2.4cm H <sub>2</sub> O	61ml	32g	200ml

<b>1941001</b>	70	<b>1841000</b>	35
Filta-Therm Plus + luer lock port		Clear-Therm + luer lock port	

## Inter-Therm HMEF<sup>®</sup>

The Inter-Therm HMEF provides both high filtration efficiency and heat and moisture performance and is provided **sterile**.

Filtration efficiency	Moisture return at: VT500ml	Resistance to flow		Compressible volume		Weight		Minimum tidal volume
		30 L/min	60 L/min	Without port	With port	Without port	With port	
99.999	32mg H <sub>2</sub> O/L	1.6cm H <sub>2</sub> O	3.0cm H <sub>2</sub> O	56ml	57ml	30g	31g	150ml

<b>1341007S</b>	125	<b>1341000S</b>	125	<b>1341580S</b>	50
Inter-Therm HMEF		Inter-Therm HMEF + luer lock port and elbow		Inter-Therm HMEF + luer lock port and Smoothbore catheter mount with double swivel elbow and double flip top cap	



## Hydro-Therm® HME

A range of dedicated HME's where filtration is not required

Moisture return at: VT 500ml HME only	Resistance at: HME only		Compressible volume HME only		Weight HME only		Minimum tidal volume HME only
	30L/min	60L/min	Without port	With port	Without port	With port	
<b>Hydro-Therm</b> 30mg H <sub>2</sub> O/L	0.7cm H <sub>2</sub> O	1.9cm H <sub>2</sub> O	15ml	16ml	11g	12g	50ml
<b>Hydro-Therm II</b> 33mg H <sub>2</sub> O/L	0.3cm H <sub>2</sub> O	1.4cm H <sub>2</sub> O	N/A	60ml	N/A	33g	200ml

1850	20	1855	20	1860	35
Hydro-Therm		Hydro-Therm + luer lock port		Hydro-Therm II + luer lock port	

## Flextube® and Smoothbore basic breathing systems

A range of basic two limb breathing systems in both Flextube and Smoothbore tubing for use with HME's or HMEF's for a passive humidification solution.

		2000
<b>2000</b>	Breathing system, 1.6m length	20
<b>2000100</b>	Anti-microbial breathing system, 1.6m length	20
<b>2200</b>	Breathing system with detachable limb, 1.6m length (for insertion of a nebuliser or MDI)	20
<b>2006</b>	Breathing system, 2.0m length	15
<b>2001</b>	Breathing system, 2.4m length	15
<b>2002</b>	Breathing system, 1.6m length with 0.8m limb	20
<b>2003</b>	Breathing system, 2.4m length with 0.8m limb	15
<b>2004</b>	Breathing system with ported Y-piece, 1.6m length with 0.8m limb	10
<b>5000</b>	<b>500000S</b> Smoothbore breathing system 1.6m length	12




## Active Humidification

If your ventilated patient requires active humidification then our range of humidification chambers combined with a choice of breathing systems in Flextube or Smoothbore are available.

### Humidification Chambers

The humidification chambers are an integral part of the breathing system and allow the system to interface with the heated humidifier base. The range consists of three chamber options which can be used with humidifier bases commonly used in intensive care units.

The chamber simply slides into position on the hot plate of the base controller allowing the inspiratory gas to pass over the heated water. These are available with breathing systems for convenience or individually if required.



#### Auto-fill humidification chamber


The auto-fill humidification chamber offers a fixed level of water within the chamber, ensuring a constant system volume. This, coupled with the strong polycarbonate body and non compressible float, ensures that adverse changes in system compliance are reduced to a minimum.

The auto-fill chamber provides optimum humidification output without compromising resistance to flow. The new dual-float, dual-valve design provides further assurance of reliability.

2310

Auto-fill humidification chamber

30




#### Manual fill humidification chamber

The manual fill humidification chamber offers a cost effective option in all areas of ventilation. This product is supplied complete with fill set and clamp in order to manually control the water level in the chamber.

2320

Manual fill humidification chamber

30



#### Low volume humidification chamber

The low volume humidification chamber is suitable for use with high frequency ventilation and many neonatal applications. The product is supplied with a fill set and clamp for manually controlling the water level in the chamber.

2330

Low volume manual fill humidification chamber

30



## Features and benefits of the auto-fill humidification chamber

### Inline filter

Prevents any debris entering the chamber.

### Protective cassette

Protects the chamber from contamination and conveniently stores the fill set before use.

### Dual valve

For added security.

### Strong polycarbonate clear material

Improved compliance characteristics. Allows for easy visual assessment of the fluid level at all times.

### Clearly visible water level indicator

Easy to see for instant accurate fluid level assessment.

### Advanced dual float design

The closed cell material of the primary float ensures a totally reliable unsinkable rigid mechanism. Whilst the secondary float provides added security.

### Shrouded heated plate cover


Prevents accidental burns when removing the chamber from the heater base.





## Flextube® heated wire breathing systems

Systems are available with single or dual heated wires and can be supplied complete with an auto-fill humidification chamber if required




<b>2025</b>	<b>Single</b> heated wire breathing system, 1.6m length with 0.8m limb	10	
<b>2025310</b>	<b>Single</b> heated wire breathing system with <b>auto-fill chamber</b> , 1.6m length and 0.8m limb	8	
<b>2026</b>	<b>Dual</b> heated wire breathing system, 1.6m length with 0.8m limb	10	
<b>2026310</b>	<b>Dual</b> heated wire breathing system with <b>auto-fill chamber</b> , 1.6m length and 0.8m limb	8	



## Smoothbore breathing systems for active humidification

Systems are available with single or dual heated wires and can be supplied complete with an auto-fill humidification chamber if required



<b>5025</b>	<b>Single</b> heated wire breathing system, 1.6m length with 0.5m limb	10	
<b>5025310</b>	<b>Single</b> heated wire breathing system <b>with auto-fill chamber</b> , 1.6m length with 0.5m limb	8	
<b>5026</b>	<b>Dual</b> heated wire breathing system, 1.6m length with 0.5m limb	10	
<b>5026310</b>	<b>Dual</b> heated wire breathing system <b>with auto-fill chamber</b> , 1.6m length with 0.5m limb	8	



## Patient connections

A wide range of patient connections are available in Flexible, Superset and Smoothbore tubing, see our product catalogue for the full range.

All port caps are retained to ensure they cannot be misplaced in use. Two varieties of flip top cap are available to allow for suctioning and the use of a fibre optic bronchoscope.

**Double Flip top' cap with seal**

First tab open - bronchoscope mode with seal

Second full cap open - suction mode

**bronchoscopy and suction ports in one patient connection**

**Flip top cap with 7.6mm port**

**or connection of all standard 7.6mm probes or for suctioning**

Double Flip top' cap with seal

22F

Covers 7.6mm and 9.5mm applications

Swivel

22M/15F swivel

**3505** Swivel elbow and double flip top cap 50

Flip top' cap with 7.6mm port

22F

170mm

Swivel

22M/15F swivel

**3516** Swivel elbow and flip top cap 50

Covers 7.6mm and 9.5mm applications

Double Flip top' cap with seal

22F

70mm compressed  
150mm extended

Swivel

22M/15F swivel

**3521** Superset with swivel elbow and double flip top cap 75

Flip top' cap with 7.6mm port

22F

70mm compressed  
150mm extended

Swivel

22M/15F swivel

**3520** Superset with swivel elbow and flip top cap 75

Double Flip top' cap with seal

22F

Swivel

22M/15F swivel

**5180** Smoothbore with swivel elbow with double flip top cap, 180mm 50

**5180000S** Sterile Smoothbore with swivel elbow with double flip top cap, 180mm 85





# Passive Humidification

For the spontaneously breathing patient with a tracheostomy or receiving oxygen therapy, both passive and active humidification options are available.

## Hydro-Trach® T Mk

The Hydro-Trach T Mk.II is a heat and moisture exchange device designed for use with spontaneously breathing patients in order to reduce loss of heat and moisture during respiration.

When a patient has a tracheostomy, the normal system of temperature and moisture maintenance is bypassed by the insertion of the tracheal tube and can lead to serious complications.


The Hydro-Trach T Mk.II has a number of unique features which make it an ideal product for prolonged use with spontaneously breathing patients - available sterile if required.



### Features and benefits

**Clear housing**  
for easy visual inspection for possible secretion build up

**Anti occlusion mechanism**  
allowing the HME element to partially dislodge in the event of total occlusion or vigorous cough



**Small and lightweight**  
reducing the pull on the patient connection

**Clipped suctioning port**  
To allow for easy suctioning without removal of the device

**An integral swivel oxygen connector**  
allowing for connection of the oxygen tube without the need of a separate oxygen adapter, which can be easily folded away when not in use

Moisture return at: VT 500ml <small>HME only</small>	Resistance at: <small>HME only</small>		Compressible volume <small>HME only</small>	Weight <small>HME only</small>	Minimum tidal volume <small>HME only</small>
	30L/min	60L/min			
26mg H <sub>2</sub> O/L	0.2cm H <sub>2</sub> O	0.7cm H <sub>2</sub> O	19ml	8g	50ml

1873	25	1874	40
1873000S - sterile	100	1874000S - sterile	30
<p>15M</p>		<p>15M O<sub>2</sub> tube not to scale</p>	
Hydro-Trach T Mk.II		Hydro-Trach T Mk.II + O <sub>2</sub> tube	



## Aquamist™ humidifier nebulisers

For the spontaneously breathing patient receiving long term oxygen therapy humidification is essential in order to bring dry oxygen gas to ambient levels of humidity.










A number of options are available. Aerosol masks, tracheostomy masks and T-pieces provide an ideal interface for the Aquamist humidifier nebuliser which has been designed to deliver accurate concentrations of humidified oxygen quietly.

<b>1509</b>	<b>1508</b>	<b>1500001</b>	1506001	2.0m Flextube
<b>1508</b>	Humidifier nebuliser without bottle	24		
<b>1509</b>	Humidifier nebuliser with bottle	25		
<b>1500</b>	Humidifier nebuliser + BOC adaptor without bottle	24		
<b>1502</b>	Humidifier nebuliser + BOC adaptor with bottle	25		
<b>1500001</b>	Humidifier nebuliser + adult aerosol mask kit without bottle + BOC adaptor + 2.0m Flextube	15		

## Tracheostomy mask and T-pieces

<b>1200</b>	<b>1984</b>	7.6mm port
<b>1200</b>	Adult tracheostomy mask	30
<b>1206</b>	Paediatric tracheostomy mask	35
<b>1980</b>	Tracheostomy T-piece	25
<b>1984</b>	Tracheostomy T-piece with 7.6mm port	25











## Aerosol masks

<b>1103</b>	<b>1108</b>	<b>1148</b>	<b>1198</b>
			
<b>1103</b>	Adult aerosol mask	70	
<b>1101</b>	Small adult aerosol mask	70	
<b>1148</b>	Paediatric aerosol mask	35	
<b>1108</b>	Adult aerosol mask with nose clip	70	
<b>1198</b>	Paediatric aerosol mask with nose clip	50	

## AquaFlow™ oxygen bubble humidifiers

Patient's receiving variable oxygen concentrations delivered via mask, or nasal cannulae can be humidified using the Intersurgical AquaFlow. This uses the bubble-through humidification process.


















The dry gas from the flowmeter is directed into the water bottle where it is broken up into small bubbles which gain humidity as they rise to the surface of the water.

<b>1521</b>	<b>1520</b>	<b>1505</b>	<b>1506</b>	<b>1507</b>
				
<b>1506</b>	Oxygen bubble humidifier without bottle	30		
<b>1505</b>	Oxygen bubble humidifier with bottle + M12 adaptor	20		
<b>1507</b>	Oxygen bubble humidifier with bottle	20		
<b>1521</b>	Homecare bubble humidifier with bottle	20		
<b>1520</b>	Homecare bubble humidifier with bottle + M12 adaptor	20		















## Oxygen masks

### Medium concentration oxygen masks

				
<b>1136</b>	<b>1146</b>		<b>1115</b>	
<b>1135</b>	Adult Eco oxygen mask with 2.1m oxygen tube		40	
<b>1136</b>	Adult Eco oxygen mask		55	
<b>1104</b>	Adult oxygen mask		70	
<b>1104001</b>	Adult oxygen mask with ear loops		70	
<b>1105</b>	Adult oxygen mask with oxygen tube		50	
<b>1106</b>	Small adult oxygen mask		70	
<b>1115</b>	Adult oxygen mask with nose clip and oxygen tube		50	
<b>1116</b>	Adult oxygen mask with nose clip		70	
<b>1140</b>	Paediatric oxygen mask		50	
<b>1146</b>	Paediatric oxygen mask with oxygen tube		50	

## Nasal cannulae

### Straight prong, curved prong, flared prong and curved flared prong options

				
Straight prong	Curved prong	Flared prong	Curved flared prong	
<b>1161</b>	Adult straight prong with tube, 1.8m length		50	
<b>1162</b>	Adult straight prong with tube, 5.0m length		20	
<b>1169</b>	Adult straight prong headset, 0.5m length		100	
<b>1165</b>	Adult curved prong with tube, 1.8m length		50	
<b>1167</b>	Adult curved/flared prong with tube, 1.8m length		50	
<b>1168</b>	Adult curved/flared prong headset, 0.5m length		100	
<b>1166</b>	Adult flared prong with tube, 1.8m length		50	
<b>1163</b>	Paediatric curved prong with tube, 1.8m length		50	
<b>1164</b>	Neonatal curved prong with tube, 1.8m length		50	