

Silente Adjustable Venturi Valve



Silente Adjustable Venturi Valve

A new design for improved performance

The new Silente incorporates all the best features from the original design while improving the performance at high flows and drastically reducing operating noise.

This new twin-jet design offers the user seven oxygen concentrations. The increased body diameter provides space for clear, widely spaced calibration and incorporates a 'pull-twist' action preventing accidental adjustment. The one-piece design eliminates the need to use more than one valve for all concentrations and there are no separate pieces to get lost or mislaid.

Ease of Use

The Intersurgical adjustable venturi valve is simple to use and has the following features and benefits:

- The appropriate flow range is selected to suit the desired concentration
- The outer ring with its pull-twist action selects the prescribed oxygen concentration
- A flow in L/min is printed above each concentration indicating the optimum flow meter setting.

Precision

The Silente delivers precise oxygen concentrations. This accuracy is independent of the oxygen flow meter setting and is achieved using just one valve.

The concentration delivered is unaffected by the patient's breathing pattern provided that the optimum flow indicated on the valve is set on the oxygen flow meter.

Features and benefits

Large, well spaced calibration settings

Easy to read, minimises errors

Pull-Twist controls

Prevent inadvertent adjustment

One piece design

No parts to change or lose

Extremely quiet

Much improved patient environment



Ordering information

Code: 1137	Box qty. 80		
Silente, Adult Adjustable Venturi Valve Kit			
Code: 1138015	Box qty. 40		
Silente, Paediatric Adjustable Venturi Valve Kit			
Code: 1198015	Box qty. 42		
Paediatric Eco mask			
Code: 1188	Box qty. 70		
Adult Eco aerosol mask			
Code: 1139	Box qty. 50		
Silente, Adjustable Venturi Valve only			