



Bronchoscopy masks

For analgesic gas administration
during bronchoscopy procedures





Bronchoscopy masks

Bronchoscopy masks for single patient use to prevent the possibility of cross infection.

The Intersurgical bronchoscopy masks are designed to allow easy insertion of a telescope probe into the nose or mouth whilst administering an analgesic gas at the same time. There are three one-way valves positioned on the mask (with the exception of the size 1 where there is only one valve), which provides easy access from either side of the patient.

The masks are principally used for bronchoscopy and endoscopy procedures, the analgesic gas is administered via the central 22mm port. These masks can also be used on patients undergoing a transesophageal endocardiology (T.O.E) procedure.


Features and benefits

- Colour-coded hook ring**
For ease of identification
- Anatomical shape**
For a precise fit.
- Distinct colour-coding**
For easy selection and integration with other Intersurgical airway products.
- 4 sizes available**
From Adult to Infant
- Soft cushion**
For patient comfort and air tight fit.
- Self-sealing valves**
To facilitate easy access.

More information?

Call **0118 9656 376**
www.intersurgical.com

Product information



Code: 1591	Box Qty. 40
Extra large adult anaesthetic face mask size 1 ■	
Code: 1592	Box Qty. 35
Large adult anaesthetic face mask size 2 □	
Code: 1593	Box Qty. 50
Medium adult anaesthetic face mask size 3 ■	
Code: 1594	Box Qty. 40
Small adult anaesthetic face mask size 4 ■	