

STARMed Ventumask

For CPAP therapy with Venturi flow driver and adjustable PEEP valve



StarMed Ventumask – CPAP therapy with Venturi flow driver and adjustable PEEP valve

StarMed Ventumask is a unique design for CPAP therapy that is simply activated by connecting to an appropriate oxygen supply.

Ventumask offers:

- Convenience
- Ease of use
- Reliability
- Immediate control of the basic therapy parameters







 
CE certified product

Table of Ventumask performance

Requirement		Setting	Flow Delivered		Requirement		Setting	Flow Delivered	
FiO ₂ (*) approx	Peep cm H ₂ O	Flow O ₂ l/min (A) + (B)	Flow approx l/min		FiO ₂ (*) approx	Peep cm H ₂ O	Flow O ₂ l/min (A) + (B)	Flow approx l/min	
40%	5	12 + 0	57		60%	5	8 + 20	60	
	7.5	12 + 0	52			7.5	8 + 20	55	
	10	12 + 0	47			10	8 + 20	50	
	12.5	12 + 0	42			12.5	9 + 16	42	
50%	5	10 + 12	58		100%	Feed with a minimum O ₂ flow of 42 l/min			
	7.5	10 + 12	54						
	10	10 + 10	47						
	12.5	10 + 10	42						

Note – The use of flow meters with an output of less than 30 litres/minute on line B, will reduce the performance of the device above 50% FiO₂

CV200M	CPAP mask with integrated Venturi flow driver and adjustable PEEP valve	Small / Medium	5	
CV200L	CPAP mask with integrated Venturi flow driver and adjustable PEEP valve	Large	5	
CV200XL	CPAP mask with integrated Venturi flow driver and adjustable PEEP valve	Extra Large	5	
STS.15/30	Dual Flow Meter STS.15/30 flow meter 15/30 L/min with O ₂ coupling		1	

Single packed. Minimum purchasing package: 5 units.

IS123_EN Issue 3 12.12



Manufactured by:
Intersurgical S.p.A. - Via T. Morandi, 12
41037 Mirandola MO - Italy
T. +39 0535 610131 - F. +39 0535 610310

UK
info@intersurgical.co.uk
Italy
info@intersurgical.it

España
info@intersurgical-es.com
Portugal
info@intersurgical.pt

Lietuva
info@intersurgical.lt
Россия
info@intersurgical.ru

France
info@intersurgical.fr
South Africa
info@intersurgical.co.za

Japan
info@intersurgical.co.jp
Nederland
info@intersurgical.nl

