

STARMed CaStar Infant

Infant hood for CPAP therapy



STARMed CaStar Infant hood for CPAP therapy

StarMed CaStar Infant


The CaStar Infant hood is a comfortable, versatile and lightweight CPAP patient interface for infants. It is designed to offer an alternative to traditional methods of CPAP delivery. The CaStar Infant hood is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in both semi-recumbent and supine positions and is available in two sizes, providing an ideal option for CPAP therapy for a range of patient sizes.

The CaStar Infant hood may be used for the treatment of a number of respiratory conditions including:



Pre and post extubation respiratory support, thoracic trauma, asthma, post operative hypoxaemia, cardiogenic and non cardiogenic pulmonary oedema.

Technical Features

- 22M/22F connectors compliant to ISO standard
- Available in 2 sizes:
 - Low internal volume 7 litres
 - High internal volume 8 litres
- Provided with built-in pressure manometer to check the presence of pressure within the hood
- Fitted with inflatable headrest
- Atraumatic harness for patient comfort
- Two sealed access ports for probes or catheters of 2.5-4.0mm Ø
- Patient access port with screw cap. Fitted with bidirectional anti-suffocation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated



Features and benefits

-   CE certified product
- Transparent**
For ease of patient monitoring
- Integrated pressure manometer**
- Patient access port with anti-suffocation valve**
- Comfortable patient harness**

Code	Description	Approx Patient weight	Box qty.
CP200INFHI/2	CaStar Infant High with anti-suffocation valve and patient access port	From 7 to 12 Kg	5
CP200INFLO/2	CaStar Infant Low with anti-suffocation valve and patient access port	From 3 to 7 Kg	5

Single packed. Minimum purchasing package: 5 units. Custom-packs are supplied on specific request of the customer.



INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS
Crane House, Molly Millars Lane,
Wokingham, Berkshire, RG41 2RZ
T. +44 (0)118 9656 300
F. +44 (0)118 9656 356

UK
info@intersurgical.co.uk
Italy
info@intersurgical.it

España
info@intersurgical-es.com
Portugal
info@intersurgical.pt

Lietuva
info@intersurgical.lt
Россия
info@intersurgical.ru

France
info@intersurgical.fr
South Africa
info@intersurgical.co.za

Japan
info@intersurgical.co.jp
Nederland
info@intersurgical.nl



IS1.20 • Issue 3_05.13

Manufactured by: Intersurgical S.p.A.