

STARMed CaStar 'R' NEXT

Hood without underarm straps for non invasive mechanical ventilation (NIV)



STARMed CaStar 'R' Next hood for non invasive ventilation (NIV)

StarMed CaStar R Next

The CaStar R Next is a comfortable, versatile and lightweight patient interface for non invasive mechanical ventilation (NIV) designed to offer an alternative to traditional methods of NIV delivery. The CaStar R Next is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in both semi-recumbent and supine positions providing an ideal option for NIV delivery.

The CaStar R Next has been designed to improve patient-ventilator interaction, by reducing the inspiratory trigger delay and increasing the rate of airway pressurisation, while avoiding the need for underarm straps.

The CaStar R Next may be used for the treatment of a number of respiratory conditions including:


Pre and post extubation respiratory support, COPD exacerbation, thoracic trauma, asthma, post operative hypoxaemia, cardiogenic and non cardiogenic pulmonary oedema.

Reference: Intensive Care Med. 2013 Spr39(4): 734-8
Bench comparative evaluation of a new generation and standard helmet for delivering non-invasive ventilation.
 Oliveri C, Costa R, Spinazzola G, Ferrone G, Longhini F, Cammarota G, Conti G, Navalesi P.

Features and Benefits



CE certified product



New fastening system

Replaces underarm straps



New sealed access

For probes or catheters

CaStar R Next hood with anti-suffocation valve, patient access port, ear plugs and measuring tape

Code	Size	Neck Circumference	Box qty.
CP238XXL/2R	Extra Extra Large	48-55cm	5
CP238XL/2R	Extra Large	42-48cm	5
CP238L/2R	Large	35-42cm	5
CP238M/2R	Medium	28-35cm	5
CP238S/2R	Small	21-28cm	5
CP238XS/2R	Extra Small	≤21cm	5

Single packed. Minimum purchasing package: 5 units. Custom-packs are supplied on specific request of the customer.

IS1.24_ROW • Issue 3_05.13



Manufactured by:
 Intersurgical S.p.A. - Via T. Morandi, 12
 41037 Mirandola MO - Italy
 T. +39 0535 610131 - F. +39 0535 610310

UK
info@intersurgical.co.uk
 Italy
info@intersurgical.it

España
info@intersurgical-es.com
 Portugal
info@intersurgical.pt

Lietuva
info@intersurgical.lt
 Россия
info@intersurgical.ru

France
info@intersurgical.fr
 South Africa
info@intersurgical.co.za

Japan
info@intersurgical.co.jp
 Nederland
info@intersurgical.nl

