


# Laryngeal Mask Airways

A comprehensive range of **single use** laryngeal mask airways for use in anaesthesia and emergency medicine



## Solus® standard laryngeal mask airway

Single use laryngeal mask airway ideal for use in anaesthesia and emergency medicine



[www.youtube.com/intersurgical](http://www.youtube.com/intersurgical)



Code	8005	8004	8003	8025	8002	8015	8001
Box qty	20	20	20	20	20	20	20
Size	5	4	3	2.5	2	1.5	1
Weight (kg)	>70	50-70	30-50	20-30	10-20	5-10	<5

## Solus® Satin laryngeal mask airway

With softer airway tube



Code	8005003	8004003	8003003
Box qty	20	20	20
Size	5	4	3
Weight (kg)	>70	50-70	30-50

## Solus® Flexible wire-reinforced laryngeal mask airway

Permits flexion without kinking



Code	8005001	8004001	8003001	8025001	8002001
Box qty	20	20	20	20	20
Size	5	4	3	2.5	2
Weight (kg)	>70	50-70	30-50	20-30	10-20

## Solus® MRI safe laryngeal mask airway

With non-ferrous valves guaranteed not to interfere with the magnet in an MRI unit



Code	8005002	8004002	8003002	8025002	8002002	8015002	8001002
Box qty	20	20	20	20	20	20	20
Size	5	4	3	2.5	2	1.5	1
Weight (kg)	>70	50-70	30-50	20-30	10-20	5-10	<5