



QuadraLite™

Anatomically-shaped with a flexible seal
in four sizes to fit all patients



Introducing QuadraLite™

Soft seal

A comfortable airtight fit with flexibility for different face sizes

Transparent shell

For excellent visibility



Colour coding

Colour matching with other Intersurgical products for easy selection

Anatomical shape

Provides a precise fit combined with lower dead space

Less means more, to the power of four

The QuadraLite is a soft and comfortable anaesthesia mask that is available in four sizes to fit all patients. The unique anatomical shape combined with its flexible seal means that only four sizes are required, unlike traditional cushion masks which often require six or seven sizes. This means you can hold less stock and still provide a complete solution for all of your patients.

The greener, safer choice

QuadraLite masks contain no PVC and are phthalate free.



PHthalate-FREE

Comfort without the cushion


QuadraLite masks do not have a cushion, but an innovative use of two different materials means the mask offers an optimum level of flexibility and comfort. A crystal-clear shell enhances patient visibility.

*QuadraLite was formerly known as EcoAnaesthesia

Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	QuadraLite masks
Neonatal	Preterm	< 2	0	 <p>Code: 7191</p>
	Term	2 – 4		
Infant	3 months	6	1	
	6 months	8		
	9 months	9		
Paediatric	1 year	10	2	
	2 years	12		
	3 years	14		
	4 years	16		
	5 years	18		
Small Adult	6 years	20	3	
	7 years	23		
	8 years	26		
	9 years	28		
	10 years	32		
	11 years	35		
Medium Adult	12 years	39	4	 <p>Code: 7193</p>
	13 years	44		
	14 years	50		
	15 years	54		
Large Adult	16 years	58	5	
	16+ years	60 – 120		
Extra Large Adult	16 + years	> 120	6	 <p>Code: 7194</p>

Ordering information



QuadraLite anaesthesia mask 			
7194	Large anaesthetic face mask with hook ring size 5-6		35
7193	Medium anaesthetic face mask with hook ring size 3-4		35
7192	Small anaesthetic face mask with hook ring size 2-3		25
7191	Infant anaesthetic face mask size 0-1		30

