



# **Sterile Smoothbore Range**

Breathing systems and catheter mounts



# Now there's more choice in our sterile range

We are pleased to announce that we have extended our range of sterile options to offer smoothbore breathing systems and catheter mounts.

Smoothbore tubing gives improved ventilatory support with it's improved flow characteristics, lower tubing compliance and clear crush-resistant tubing.

Other benefits include:

- Minimal water collection within the system and improved drainage
- Soft flexible end connectors to overcome non-standard fittings
- Lower system volume

**Even more choice**  
If you require a breathing system that is not shown here, please contact your Customer Service contact for details of alternative options.

## Ordering information

<b>512000S</b>	Smoothbore catheter mount with swivel elbow and double flip top cap, 120mm	50
<b>518000S</b>	Smoothbore catheter mount with swivel elbow and double flip top cap, 180mm	50
<b>525000S</b>	Smoothbore catheter mount with swivel elbow and double flip top cap, 250mm	50

<b>500900S</b>		
<b>500000S</b>	Basic breathing system for passive humidification, 1.6m length	12
<b>500700S</b>	Breathing system with ported Y-piece and water trap, 1.6m length	10
<b>500900S</b>	Breathing system with ported Y-piece and two water traps, 1.6m length with 0.5m limb	10