

# nFlow<sup>™</sup>

Infant nasal CPAP system



*Designed to care*

# nFlow™

## Designed to care

The smallest patients often need the most care. Safety and comfort are vital to help infants with their earliest steps in life.

## A new choice for infant CPAP

The new Intersurgical nFlow system offers a complete care solution for infant CPAP with established 'gas flip' technology, providing a low 'work of breathing' option as an alternative to ventilator CPAP.



## The nFlow system


The nFlow system (code 4700) is supplied assembled with generator, heated wire limb, 10mm adjustable Superset style exhaust tube for efficient gas removal, and phthalate free fresh gas and pressure monitoring lines.

## A choice of patient interfaces and bonnets

Patient interfaces and bonnets are supplied separately to avoid waste. Use the size guide provided to choose the correct option for your patient.


## Designed for compatibility

nFlow is compatible with all the common flow drivers and humidifiers commonly found in the neonatal intensive care areas.



## Choice made easy

The nFlow system includes two sizing guides making it really simple to choose the patient interface and bonnet.



**Patient interface**  
Masks and prongs are available separately in three size options. Made from soft, comfortable silicone to ensure optimal patient fit with minimal wastage.

**12 bonnets**  
Available separately in twelve colour-coded sizes. Each bonnet includes three Velcro fasteners and secure attachment points.

### Masks



### Prongs




### Bonnets



# How to order nFlow™

We've made ordering easy. To make a complete system you can choose one product code from each group marked **A to D**.

At any time you can order separate items according to your needs.



## A nFlow system



<b>4700000</b>	nFlow infant CPAP system with heated wire 1.6m length	<b>12</b>
----------------	---	-----------

## B Bonnets



<b>4707009</b>	Size 9 (40 - 42 cm)	Navy		<b>10</b>
<b>4707008</b>	Size 8 (38 - 40 cm)	Turquoise		<b>10</b>
<b>4707007</b>	Size 7 (36 - 38 cm)	Orange		<b>10</b>
<b>4707006</b>	Size 6 (34 - 36 cm)	Red		<b>10</b>
<b>4707005</b>	Size 5 (32 - 34 cm)	Green		<b>20</b>
<b>4707004</b>	Size 4 (30 - 32 cm)	Light Orange		<b>20</b>
<b>4707003</b>	Size 3 (28 - 30 cm)	Light Blue		<b>20</b>
<b>4707002</b>	Size 2 (26 - 28 cm)	Yellow		<b>20</b>
<b>4707001</b>	Size 1 (24 - 26 cm)	Brown		<b>20</b>
<b>4707110</b>	Size 0 (22 - 24 cm)	Pink		<b>20</b>
<b>4707100</b>	Size 00 (20 - 22 cm)	Grey		<b>20</b>
<b>4707000</b>	Size 000 (18 - 20 cm)	White		<b>20</b>

## C Nasal prongs



<b>4703000</b>	Large	Purple		<b>30</b>
<b>4702000</b>	Medium	Blue		<b>30</b>
<b>4701000</b>	Small	Pink		<b>30</b>

## D Nasal masks



<b>4706000</b>	Large	Purple		<b>30</b>
<b>4705000</b>	Medium	Blue		<b>30</b>
<b>4704000</b>	Small	Pink		<b>30</b>

IS1.16 • Issue 3 05.12