



STR/017
CPAP Nebuliser Accessory, complete with mouthpiece to provide an aerosol medication without removing the CaStar hood from the patient's head.



STR/014
22F/22F Straight Adaptor can be used with the CaStar R to enable a PEEP valve to be fitted to the exhalation port so that it can be used as a CPAP Hood.



SPV12F
22F Cap with leak. For use with the CaStar R (depending upon the type of breathing circuit used)



STR/019
Cap. For use with the CaStar R (depending upon the type of breathing circuit used)



CP903
Cuff Inflation Bulb for CaStar R. For convenient inflation of the cuff within the CaStar R



CP900
Underarm Strap.