



**STARMed**  
medical disposables

# Ventukit Hood

for CPAP therapy



# Ventukit for CPAP therapy

### The Ventukit contains:

A CaStar CPAP hood, Venturimeter and silencer. CPAP hood with patient inspection port. Bi-directional anti-suffocation valve. Two sealed ports for probes or catheters of 3.5-7.0mm Ø. Biocompatible and latex-free.



Integrated Venturimeter for FiO<sub>2</sub> from 40% to 100%, equipped with silencer to reduce the noise inside the hood by up to 10 dB, and a port for instantaneous FiO<sub>2</sub> control.



Integrated pressure gauge and adjustable PEEP valve.



'No-touch' earplugs to further reduce the noise during extended therapies.



The Ventukit circuit is connected to a dual flowmeter connected to the O<sub>2</sub> wall outlet (probe supplied separately).

The therapy is activated simply by connecting The Ventukit to an oxygen source commonly available.

The Ventukit offers:

- Convenience
- Ease of use
- Reliability
- Immediate control of the basic therapy parameters.

FiO <sub>2</sub> % <i>approx</i>	O <sub>2</sub> L/min <b>(A) + (B)</b>	PEEP cm H <sub>2</sub> O	Flow L/min <i>approx</i>
40	12 + 0	7.5	48
50	10 + 7	7.5	47
60	9 + 14	7.5	48

**Table of device performance with a PEEP of 7.5cm H<sub>2</sub>O.**

For other PEEP settings see the instructions for use provided with the product.

Model	Description	Neck Size	Quantity	Icon
<b>CV100S/2</b>	Small Ventukit CPAP Hood with safety valve and patient access port	<b>27 - 34cm neck</b>	5	
<b>CV100M/2</b>	Medium Ventukit CPAP Hood with safety valve and patient access port	<b>34 - 41cm neck</b>	5	
<b>CV100L/2</b>	Large Ventukit CPAP Hood with safety valve and patient access port	<b>40 - 47cm neck</b>	5	
<b>CV100XL/2</b>	XL Ventukit CPAP Hood with safety valve and patient access port	<b>≥45cm neck</b>	5	