


# Packaging Symbols Explained

Symbol	Meaning	Symbol	Meaning
	Catalogue number		Sterilized using Ethylene Oxide
	Batch Code		Do not re-sterilize
	Use by date		Do not use if package is damaged
	Caution, consult accompanying documents		Non-sterile
	Do not re-use		Upper limit of temperature
	Date of manufacture		Temperature limitation
	Manufacturer		MRI compatible – will not interfere with the magnet in an MRI scanner
	CE mark. Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC.		Irritant – may cause inflammation to the skin or other mucous membranes
	CE mark including SGS (notified body) identification number. Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC.		Keep away from sunlight
	Product does not contain latex		Keep away from sunlight
	Product does not contain latex		This way up
	Product does not contain PVC		Recyclable
			Do not open packaging using a knife