

August 2009

Intersurgical Breathing Filters

In the light of the current Swine Flu (H1N1) global pandemic, the question has been raised concerning the suitability of Intersurgical breathing filters at preventing the passage of influenza viruses in respiratory systems.

The range of Intersurgical breathing filters has been independently tested to be highly efficient at preventing the passage of a variety of viral species of varying particle sizes greater than 23nm¹ including *MS-2 coliphage*¹, *Phi-x174 bacteriophage*² and *Hepatitis C*³.

Influenza viruses, including the H1N1 strain, range in size from 80-120nm⁴. On this basis Intersurgical filters should also be highly efficient at preventing the passage of the Swine Flu virus.

Glyn Jones

Filter Product Specialist
Intersurgical

Ref

1. DUBOVI, E. J. and AKERS, T. G. (1970). Airborne stability of tailless bacterial viruses S-13 and MS-2. *Appl. Microbiol.* 19, 624-628.
2. Nelson Laboratories inc, Salt Lake City, Utah, USA SOP/ARO/018G.1
3. YUASA, T., ISHIKAWA, G., MANABE, S., SEKIGUCHI, S., TAKEUCHI, K, and MIYAMURA, T. (1991), The particle size of hepatitis C virus estimated by filtration through microporous regenerated cellulose fibre. *J. Gen. Virol.* 72, 2021-2024.
4. International Committee on Taxonomy of Viruses. ICTVdB Management (2006). 00.046.0.01. Influenzavirus A. In: *ICTVdB - The Universal Virus Database*, version 4. Büchen-Osmond, C. (Ed), Columbia University, New York, USA