

## Nasal Cannulae

A range of cannulae for adults, paediatrics and neonates



# Nasal Cannulae

Intersurgical is pleased to offer a complete range of nasal cannulae. This comprehensive range is entirely latex-free and includes designs suitable for use with adults, paediatrics and neonates.

## Choice of options

Product choices are available to suit neonatal, paediatric and adult patients. The following prong shapes are also available:

- Curved
- Flared
- Curved/Flared
- Straight


All oxygen tubing in the range is manufactured in a 'multi-channel' design, which permits oxygen flow even if the tubing is 'kinked'.

## Now available with ear guards


Ear guards are an additional comfort benefit that come pre-attached to the cannulae designed for patient satisfaction.

## Also available in Satin


Cannulae made from Intersurgical's soft satin material, proven in testing to create 20% less friction<sup>1</sup> to the patients face and ears compared to a standard cannulae




**Curved prongs:**  
Improved anatomical fit



**Flared prongs:**  
Slows down the flow of oxygen



**Curved / Flared prongs:**  
Combines the benefits of both the above designs



**Straight prongs**

## References

1. Annan. J. 69381 15.04.13

## Features and Benefits

### Multi-Channel oxygen tubing

Allows oxygen to continue to flow if the tubing becomes kinked

### Latex free

Suitable for patients with latex allergy

### Over-ear design

Secure positioning  
(Ear guards available for added comfort)

### Non-slip toggle

Enables user to adjust for best possible fit

### Clear material

Unobtrusive  
(Satin material available also)

### Soft prongs

Minimise irritation

### Location saddle

Assists in the optimum positioning of the prongs

### Adult headset version


Ideal for domiciliary use



## Ordering information

## Adult nasal cannulae

1161000



|                                   |                                                                     |     |              |
|-----------------------------------|---------------------------------------------------------------------|-----|--------------|
| 1161000                           | Adult, nasal cannulae with straight prongs and tube, 1.8m           |     | Box qty. 50  |
| 1161004                           | Adult, nasal cannulae with straight prongs and tube, 2.1m           |     | Box qty. 50  |
| 1162000                           | Adult, nasal cannulae with straight prongs and tube, 5m             |     | Box qty. 20  |
| 1169000                           | Adult, nasal cannulae with straight prong and headset, 0.5m         |     | Box qty. 100 |
| 1161020                           | Adult, Satin nasal cannulae with straight prongs and tube, 1.8m     | NEW | Box qty. 50  |
| 1161021                           | Adult, Satin nasal cannulae with straight prongs and tube, 2.1m     | NEW | Box qty. 50  |
| <div>1165002</div> <div>NEW</div> |                                                                     |     |              |
| 1165002                           | Adult, nasal cannulae with curved prongs, ear guards and tube, 2.1m | NEW | Box qty. 50  |
| 1165003                           | Adult, nasal cannulae with curved prongs, ear guards and tube, 5m   | NEW | Box qty. 40  |
| 1165000                           | Adult, nasal cannulae with curved prongs and tube, 1.8m             |     | Box qty. 50  |
| 1165020                           | Adult, Satin nasal cannulae with curved prongs and tube, 1.8m       | NEW | Box qty. 50  |
| 1165021                           | Adult, Satin nasal cannulae with curved prongs and tube, 2.1m       | NEW | Box qty. 50  |
| 1166000                           | Adult, nasal cannulae with flared prongs and tube, 1.8m             |     | Box qty. 50  |
| 1167000                           | Adult, nasal cannulae with curved/flared prongs and tube, 1.8m      |     | Box qty. 50  |
| 1168000                           | Adult, nasal cannulae with curved/flared prongs and headset, 0.5m   |     | Box qty. 100 |

## Neonatal and paediatric nasal cannulae

|         |                                                        |             |
|---------|--------------------------------------------------------|-------------|
| 1164000 | Neonatal, nasal cannulae curved prong with 2.1m tube   | Box qty. 50 |
| 1163000 | Paediatric, nasal cannulae curved prong with 2.1m tube | Box qty. 50 |

## Phalate free range of nasal cannulae



|         |                                                             |             |
|---------|-------------------------------------------------------------|-------------|
| 1160000 | Premature, nasal cannulae with curved prongs and tube, 2.1m | Box qty. 50 |
| 1160001 | Neonatal, nasal cannulae with curved prongs and tube, 2.1m  | Box qty. 50 |
| 1160002 | Infant, nasal cannulae with curved prongs and tube, 2.1m    | Box qty. 50 |

IS10.8 NL • Issue 6 09.13